

Cofresh

Polyethylene Glycol 400, Glycerine & Hypromellose

Composition

Cofresh Sterile Eye Drops: Each ml contains Polyethylene Glycol (PEG 400) BP 10mg, Glycerine BP 2 mg & Hypromellose BP 3.6 mg.

Description

Cofresh is an almost colorless, clear aqueous lubricant sterile ocular solution administered to the eye(s) in the form of drops. After topical application of this eye drop, it spreads rapidly over the conjunctiva and cornea forming a lubricating and protective film. Polyethylene Glycol 400 increases the stability of tear film and due to greater muco-adhesive properties hypromellose helps the tear film to attach with the ocular surface. Glycerin restores moisture to the dry eye. Electrolytes which are similar to the natural tears maintain good corneal epithelial surface.

Indications

It is indicated for the symptomatic relief of burning, irritation, and discomfort due to dryness of the eye or exposure to wind or sun. It is also indicated for the treatment of keratoconjunctivitis sicca.

Dosage & Administration

Instill 1 drop in the affected eye(s) 3 times daily or as needed.

Side Effect

Occasionally mild, transient burning or sticky sensation and very rarely irritation or hypersensitivity reactions reported. Blurred vision after application may occur.

Contraindication

This eye drop is contraindicated in patients with known hypersensitivity to any ingredient of the product.

Precaution

Remove contact lenses before using this drug. Do not use if this solution changes color or becomes cloudy. Don't touch tip of container to any surface to avoid contamination and replace cap after each use.

Use in Pregnancy and Lactation

There is no experience regarding the safety of this eye drops in human pregnancy or lactation. If pregnant or breast-feeding, ask a doctor before use.

Drug interaction

If several medicines are to be administered to the eye, there should be an interval of at least 5 minutes between each application.

Storage

Store at room temperature. It is desirable that the contents should not be used more than 4 weeks after first opening of the bottle.

Presentation

Cofresh Sterile Eye Drops: Each plastic dropper bottle contains 10 ml of sterile lubricant eye drops.



Manufactured by
The **IBN SINA** Pharmaceutical Industry Ltd.
Shafipur, Kaliakoir, Gazipur, Bangladesh.